

HOW TO ADD RACE DENTAL TO YOUR MEDIT SCANNER



Thank you for choosing Race Dental for laboratory needs.

Please find below instructions on how to add Race Dental to your scanner:

- 1 Log in to your Medit account
<https://www.meditlink.com/login>



The screenshot shows the MEDIT Link login interface. It features a blue header with the text "MEDIT Link". Below the header are input fields for "ID" (with "Email Address" written below it), "Password", and another "Password" field with an eye icon. There is a checkbox for "Keep me logged in." and a link for "Need help?". A prominent blue "Log In" button is centered below the fields. At the bottom, there is a link "Don't have an account? Sign up." and a language dropdown menu set to "English (United States)".

- 2 On the left side panel, select **'Partners'** and search for **'Race Dental'**. Click on **'Request Partnership'**.

A request will be sent to Race Dental's Digital Team to approve your request.

