

## Premium Restorations

	 <b>Crystalite™</b>	 <b>e.max® CAD</b>	 <b>Pressed e.max®</b>	 <b>Zirconia (Layered)</b>
Material	Monolithic Super Translucent Zirconia	Ivoclar™ e.max® lithium disilicate	Ivoclar™ e.max® lithium disilicate	Zirconia (layered) coping with porcelain overlay
Method of manufacture	Digitally designed and manufactured	Digitally designed and manufactured	Waxed & pressed by hand	Digitally designed & manufactured (Framework) & built up by hand
Aesthetics	★★★★★	★★★★★	★★★★★	★★★★★
Shades / colouring	No internal effects standard shades	Customised shading	Customised shading	Internal effects customised shading multi-layered
Strength	★★★★★	★★★★★	★★★★★	★★★★★
Add-on options <small>(fit to denture, attachments etc. Charges apply)</small>	Fit to denture	✗	✗	Fit to denture, rests and ledges
Benefits	High durability extreme strength metal-free, exceptional accuracy	Premium aesthetics excellent fit metal-free	Premium aesthetics metal-free	Good fit metal free
Recommended use	Full mouth, single crowns	Full mouth, crowns, inlays/onlays, veneers, 3/4 crowns	Full mouth, crowns, inlays/onlays, veneers, 3/4 crowns	Full mouth, crowns & bridges
Warranty	15 years	5 years	1 year	8 years

**We recommend Digital Impressions!**  
All intraoral scanners accepted







Fast 5 lab days turnaround



100% Locally Made

## Standard Restorations

	 <b>Opalite™</b>	 <b>Full Gold Crown</b>	 <b>Enamic</b>	 <b>PFM</b>
Material	Monolithic Translucent Zirconia	Non-precious (gold-coloured) or precious (gold) alloy CoCr (silver coloured)	Hybrid Ceramic	Non-precious or precious alloy coping with porcelain overlay
Method of manufacture	Digitally designed and manufactured	Traditionally-cast	Digitally designed and manufactured	Traditionally-cast & built up
Aesthetics	★★★★★	★☆☆☆☆	★★★★★	★★★★★
Shades / colouring	No internal effects standard shades	Not required	No internal effects standard shades	No internal effects standard shades
Strength	★★★★★	★★★★★	★★★★★	★★★★★
Add-on options <small>(fit to denture, attachments etc. Charges apply)</small>	Fit to denture	✓	✗	✓
Benefits	High durability extreme strength metal-free, exceptional accuracy	Exceptional strength	Minimal invasive, shock absorbing	Traditional method choice of alloys
Recommended use	Full mouth, crowns & bridges	Posterior teeth	Full mouth, crowns, inlays/onlays, implants	Full mouth, crowns & bridges
Warranty	15 years	6 years	5 years	6 years

**We recommend Digital Impressions!**  
All intraoral scanners accepted



Fast 5 lab days turnaround



100% Locally Made